


## How does it work?



### Ab & pAG-MNase binding:

1. Cells are immobilized on Concanavalin A Magnetic Beads to allow for subsequent buffer and reagent exchanges.
2. Cell membranes are then permeabilized with digitonin to facilitate the entry of primary antibody and pAG-MNase fusion enzyme into the cell nuclei.
3. The target-specific primary antibody recruits the pAG-MNase to the chromatin through protein-protein interactions between the antibody and the pAG domain of the fusion enzyme.

### MNase Digestion:

4. The addition of  $Ca^{2+}$  activates the pAG-MNase, which gently cleaves and liberates the desired chromatin fragments, allowing them to diffuse away from the genomic chromatin, out of the cell, and into the supernatant.

### DNA Purification:

5. DNA is purified using DNA purification spin columns or phenol/chloroform extraction followed by ethanol precipitation. The purified, enriched DNA is then identified and quantified using qPCR or NG-seq.

**Try it now and get 25% off on the following products:**

- CUT&RUN Assay Kit #86652
- CUT&RUN pAG-MNase and Spike-In DNA #40366