

ECB0020L

Leibovitz's L-15 medium w/o L-glutamine

Description	Leibovitz's L-15 medium w/o L-glutamine
Format	500 ml
Producer	Euroclone
Code Producer	

Applications

Suitable for established cell lines, for example Hep-2, monkey kidney cells
Qualified for primary explants from embryonic and adult human tissue
Cultivation of many viruses
Cultivation of neuronal cells

Features

Contains no sodium hydrogen carbonate
Supports the cell growth in CO₂-free systems
Buffered by the supplementation of salts, free basic amino acids and galactose

Product Specifications:

CO ₂ -Concentration, optimum	0.0 %
pH	7.0 – 7.5
Osmolality	280 – 340 mOsmol/kg
Endotoxin	≤ 1 EU/ml
Cell Culture	tested
Sterility	tested
Storage	+2°C to +8°C
Stability	18 months

Euroclone S.p.A.

