

ECB4053L

Dulbecco's Phosphate Buffered Saline

Description	Dulbecco's Phosphate Buffered Saline
Format	500 ml
Producer	Euroclone
Code Producer	

General Features:

DPBS with calcium and Magnesium

pH : 7.3 ± 0.3

Osmolality: 289 mOsm/kg $\pm 10\%$

Colour : Clear solution

Endotoxin : < 1 EU/ml

Storage conditions : $> 10^\circ\text{C}$ in the dark

Shelf life : 48 months

Sterility tests :

- Bacteria in aerobic and anaerobic conditions
- Fungi and yeasts

Recommended use :

- Respect storage conditions of the product
- Do not use the product after its expiry date
- Store product in an area protected from light
- Manipulate the product in aseptic conditions (under laminar air flow)
- Wear clothes adapted to the manipulation of the product to avoid contamination

The product is intended to be used in vitro, in laboratory only. Do not use it in therapy, human or veterinary applications.

Euroclone S.p.A.





BioNordika Sweden AB
Berzelius Väg 13, 171 65 Solna, Sweden
Tel: +46 (0)8 30 60 10 info@bionordika.se www.bionordika.se