

ECB7501L

Dulbecco's Modified Eagle's Medium High Glucose with sodium pyruvate, without L-Glutamine

Description	Dulbecco's Modified Eagle's Medium High Glucose with sodium pyruvate, without L-Glutamine
Format	500 ml
Producer	Euroclone
Code Producer	

General Features:

Liquid medium for cultivation of adherent cells and / or in suspension, with 4500 mg/l of glucose

pH : 7.3 ± 0.3
Osmolality: 330 mOsm/kg ±10 %
Colour :Red solution
Endotoxin : < 1 EU/ml
Storage conditions : +2°C to +8°C in the dark
Shelf life : 24 months

Sterility tests :

- Bacteria in aerobic and anaerobic conditions
- Fungi and yeasts

Recommended use :

- Respect storage conditions of the product
- Do not use the product after its expiry date
- Store product in an area protected from light
- Manipulate the product in aseptic conditions (under laminar air flow)
- Wear clothes adapted to the manipulation of the product to avoid contamination

The product is intended to be used in vitro, in laboratory only. Do not use it in therapy, human or veterinary applications.



Euroclone S.p.A.

Francesco Rusconi



BioNordika Sweden AB
Berzelius Väg 13, 171 65 Solna, Sweden
Tel: +46 (0)8 30 60 10 info@bionordika.se www.bionordika.se