

## ECB7504L

Dulbecco's Modified Eagle's Medium High Glucose, with Sodium pyruvate, w/o L-Glutamine and Phenol Red

<b>Description</b>	Dulbecco's Modified Eagle's Medium High Glucose, with Sodium pyruvate, w/o L-Glutamine and Phenol Red
<b>Format</b>	500 ml
<b>Producer</b>	Euroclone
<b>Code Producer</b>	

### General Features:

Liquid Dulbecco's Modified Eagle's Medium High Glucose, with Sodium pyruvate, w/o L-Glutamine and Phenol Red

pH : as reported in CoA  
Osmolality: as reported in CoA  
Colour :clear

Storage conditions : +2°C to +8°C in the dark  
Shelf life : 24 months

Sterility tests :  
-Bacteria and Fungi

### Recommended use :

- Respect storage conditions of the product
- Do not use the product after its expiry date
- Store product in an area protected from light
- Manipulate the product in aseptic conditions (under laminar air flow)
- Wear clothes adapted to the manipulation of the product to avoid contamination

The product is intended to be used in vitro, in laboratory only. Do not use it in therapy, human or veterinary applications.



Euroclone S.p.A.

*Franco Rusconi*



BioNordika Sweden AB  
Berzelius Väg 13, 171 65 Solna, Sweden  
Tel: +46 (0)8 30 60 10 info@bionordika.se www.bionordika.se