

ECB9006L

RPMI 1640 w / o L-Glutamine Sterile Filtered

Description	RPMI 1640 w / o L-Glutamine Sterile Filtered
Format	500 ml
Producer	Euroclone
Code Producer	

General Features:

RPMI 1640 Liquid medium without L-Glutamine and with Phenol Red; for cultivation of adherent cells and/or in suspension.

pH : 7.3 ± 0.3

Osmolality: 278 mOsm/kg $\pm 10\%$

Colour : Red, Orange solution

Endotoxin : < 1 EU/ml

Storage conditions : $+2^{\circ}\text{C}$ to $+8^{\circ}\text{C}$ in the dark

Shelf life : 24 months

Sterility tests :

- Bacteria in aerobic and anaerobic conditions
- Fungi and yeasts

Cell growth test :

Medium tested for the ability to support SP2/0-Ag14 cell growth

Recommended use :

- Respect storage conditions of the product
- Do not use the product after its expiry date
- Store product in an area protected from light
- Manipulate the product in aseptic conditions (under laminar air flow)
- Wear clothes adapted to the manipulation of the product to avoid contamination

The product is intended to be used in vitro, in laboratory only. Do not use it in therapy, human



or veterinary applications.

Euroclone S.p.A.

A handwritten signature in black ink, appearing to read "Franco Rusconi".