

ECM0011B

Gentamicin 10 mg/ml sulphate

Description Gentamicin 10 mg/ml sulphate

Format 10 ml

Producer Euroclone

Code Producer

General features:

Antibiotics, combined with good sterile technique, help to prevent microbiological contamination, Gram-negative and Gram-positive bacteria, mycoplasma.

Recommended use :

Use in cell culture applications at 5 ml/l. This concentration is for tissue culture media containing serum ; serum-free media generally require lower concentration.

Sterility tests :

- Bacteria in aerobic and anaerobic conditions
- Fungi and yeasts

Storage conditions : -20°C

Shelf life : 24 months

Euroclone S.p.A.

