



ECM0022L

BME with Earle's Salts

Description BME with Earle's Salts

Format 500 ml

Producer Euroclone

Code Producer

Applications:

Medium suitable for cell cultures even of immortalized cells.

Specifications:

Advised for HeLa cell cultures.

pH: 7.0-7.6

Osmol : 274-294 mOsmol/Kg

Endotoxins: < 1 EU/ml

Tested for aerobic and anaerobic bacteria, fungi and yeasts.

Conservation: +2C°-+8C°

Stability: 24 months

Euroclone S.p.A.

A handwritten signature in black ink, appearing to read "Franco Rusconi".