

ECM0090L

Dulbecco's Modified Eagle's Medium/ Ham's F-12 without L-Glutamine

Description	Dulbecco's Modified Eagle's Medium/ Ham's F-12 without L-Glutamine
Format	500 ml
Producer	Euroclone
Code Producer	

General Features:

Liquid medium for growing cells, w/o L-Glutamine

Theoretical pH : 7.4 + 0.3

Osmolality : 310 mOsm/kg + 10 %

Colour : Red solution

Storage conditions : +2°C to +8°C in the dark

Shelf life : 24 months

Sterility tests :

- Bacteria in aerobic and anaerobic conditions
- Fungi and yeasts

Endotoxin : <1 EU/ml

Cell growth test :

Medium tested for the ability to support cell growth with Hela cell line.

Preservation of Product:

The medium must be stored at 2 ° - 8 ° C away from light. Under these conditions the product may be maintained for 24 months from date of manufacture.

Composition : Displayed on web site

The product is intended to be used in vitro, in laboratory only. Do not use it in therapy, human or veterinary applications.



Uses :
Supplements, such as antibiotics, should be added as sterile supplements to the medium.
Add 12.5 ml/l of L-Glutamine 100X, 200mM before using this medium.

Euroclone S.p.A.