

## ECM0095L

Dulbecco's Modified Eagle's Medium F-12 with L-Glutamine with 25 mM Hepes

<b>Description</b>	Dulbecco's Modified Eagle's Medium F-12 with L-Glutamine with 25 mM Hepes
<b>Format</b>	500 ml
<b>Producer</b>	Euroclone
<b>Code Producer</b>	

### General features:

Liquid medium (1X) for cell cultures containing, among other components, L-Glutamine and 25mmol / liter of Hepes

### Quality Control tests:

- Sterility test of the solution by microbiological screening in appropriate media for bacteria, such as thioglycollate broth and tryptone-soy broth.
- Levels of endotoxin by LAL test.
- pH test with potentiometric measurement at room temperature
- Osmolality measurement through cryoscopic lowering analysis
- Buffer capacity test by titration with a standard acid solution.
- Cell growth promotion test through cultures of cells

Formulation available on request.

Storage conditions: the medium must be stored at +2/+8°C, away from light.

The product is intended to be used in vitro, in laboratory only. Do not use it in therapy, human or veterinary applications.

Euroclone S.p.A.





BioNordika Sweden AB  
Berzelius Väg 13, 171 65 Solna, Sweden  
Tel: +46 (0)8 30 60 10 info@bionordika.se www.bionordika.se