



## ECM0101L

# Dulbecco's Modified Eagle's Medium High Glucose w/o Sodium Pyruvate w/o L- Glutamine

<b>Description</b>	Dulbecco's Modified Eagle's Medium High Glucose w/o Sodium Pyruvate w/o L-Glutamine
<b>Format</b>	500 ml
<b>Producer</b>	Euroclone
<b>Code Producer</b>	

Liquid medium (1X) for cell cultures containing, among other components, 4500 mg/l di D-Glucose and 110 mg/l, without sodium pyruvate-

Quality control:

Microbiological screening in specific bacteria and mycoplasma mediums.

LAL tested for endotoxins

pH checked with potenziometric measurement at room temperature and osmolality with cryoscopic analysis

Tested with L929 cell cultures.

Conservation:

2° - 8°C

24 from production date

Euroclone S.p.A.