



ECM0135L

Ham's Nutrient Mixture F-12 with L-Glutamine

Description	Ham's Nutrient Mixture F-12 with L-Glutamine
Format	500 ml
Producer	Euroclone
Code Producer	

Liquid medium (1X) for cell cultures containing, among other components, L-Glutamine and Sodium Bicarbonate.

Quality control:

Microbiological screening in specific bacteria and mycoplasma mediums.

LAL tested for endotoxins

pH checked with potenziometric measurement at room temperature and osmolality with cryoscopic analysis

Tested with Vero and G7H9 cell cultures.

Conservation:

2° - 8°C

24 from production date

Euroclone S.p.A.