

## ECM0210L

McCoy's 5A w/ L-Glutamine

<b>Description</b>	McCoy's 5A w/ L-Glutamine
<b>Format</b>	500 ml
<b>Producer</b>	Euroclone
<b>Code Producer</b>	

### General Features:

Liquid medium for cell cultures, containing L-Glutamine.

### Quality Control tests:

- Sterility test of the solution by microbiological screening in appropriate media for bacteria, fungi and yeasts.
- Levels of endotoxin by LAL test
- pH test with potentiometric measurement at room temperature
- Osmolality measurement through cryoscopic lowering analysis
- Buffer capacity test by titration with a standard acid solution.
- Cell growth promotion test through cultures of cells.

### Storage conditions :

The medium must be stored at + 4 ° C away from light sources. Do not use the product after its expiry date

Euroclone S.p.A.

