

ECM0320L

Dulbecco's Modified Eagle's Medium Low Glucose with L-Glutamine with NaPyruvate

Description	Dulbecco's Modified Eagle's Medium Low Glucose with L-Glutamine with NaPyruvate
Format	500 ml
Producer	Euroclone
Code Producer	

Liquid medium (1X) for cell cultures suitable for in suspension or adhesive cells.

Quality control:

Microbiological screening in specific bacteria and mycoplasma mediums.

LAL tested for endotoxins

pH checked with potenziometric measurement at room temperature and osmolality with cryoscopic analysis.

Conservation:

2° - 8°C

12 months from production date

Euroclone S.p.A.

