

# ECM0430L

## MEM with Earle's Salts, w/o L-Glutamine, with NEAA

<b>Description</b>	MEM with Earle's Salts, w/o L-Glutamine, with NEAA
<b>Format</b>	500 ml
<b>Producer</b>	Euroclone
<b>Code Producer</b>	

### General Features:

Liquid medium MEM with Earle's Salts, w/o L-Glutamine, with NEAA, for cultivation a variety of normal and transformed cells in culture.

pH :  $7.3 \pm 0.3$

Osmolality: 290 mOsm/kg  $\pm 10\%$

Colour :Red/Orange solution

Endotoxin : < 1 EU/ml

Storage conditions : +2°C to +8°C in the dark

Shelf life : 24 months

### Sterility tests :

- Bacteria in aerobic and anaerobic conditions
- Fungi and yeasts

### Recommended use :

- Respect storage conditions of the product
- Do not use the product after its expiry date
- Store product in an area protected from light
- Manipulate the product in aseptic conditions (under laminar air flow)
- Wear clothes adapted to the manipulation of the product to avoid contamination

The product is intended to be used in vitro, in laboratory only. Do not use it in therapy, human



or veterinary applications.

Euroclone S.p.A.

A handwritten signature in black ink, appearing to read "Franco Rossetti".