

ECM0505L

RPMI 1640 w/o L- Glutamine w/o Phenol red

Description	RPMI 1640 w/o L- Glutamine w/o Phenol red
Format	500 ml
Producer	Euroclone
Code Producer	

Liquid medium (1X) for cell cultures suitable for both in suspension and adhesive cell lines.

Quality control:

Microbiological screening in specific bacteria and mycoplasma mediums.

LAL tested for endotoxins

pH checked with potenziometric measurement at room temperature and osmolality with cryoscopic analysis

Tested with murine myeloma cell lines.

Conservation:

2° - 8°C

24 months from production date

Euroclone S.p.A.

