

ECM0512L

Hank's Balanced Salts Solution w Ca/Mg w Na Bicarbonate w/o Phenol Red

Description Hank's Balanced Salts Solution w Ca/Mg w Na Bicarbonate w/o Phenol Red

Format 500 ml

Producer Euroclone

Code Producer

Balanced saline solution 1x containing::8000 mg/L NaCl 400 mg/L KCl, 60 mg/L KH₂PO₄, 1000 mg/L glucose, 47,88 mg/L Na₂HPO₄, 350 mg/L di sodium bicarbonate, 97,67 mg/L magnesium sulfate, 185 mg/L calcium chloride

Quality control:

Microbiological screening in specific bacteria and mycoplasma mediums.

LAL tested for endotoxins

pH checked with potenziometric measurement at room temperature and osmolality with cryoscopic analysis

Tested with murine NHF cell lines.

Conservation:

RT for 24 months from production date

Euroclone S.p.A.

