



## ECM0542D

Sodium Pyruvate 100 mM

**Description** Sodium Pyruvate 100 mM

**Format** 100 ml

**Producer** Euroclone

**Code Producer**

Saline solution 1X for cell cultures, (11g/L sodium pyruvate).  
Sterilized by filtration.

Quality control:

Microbiological screening in specific bacteria and mycoplasma mediums.

LAL tested for endotoxins

pH checked with potenziometric measurement at room temperature and osmolality with cryoscopic analysis

pH: pH 7 +/- 1

Osmol 200 +/- 10% mOsmol/kg

Endotoxins <1 EU/ml

Conservation:

4°C for 24 months from production date

Euroclone

Euroclone S.p.A.