

ECM0556D

MEM Vitamins 100X w/o L-Glutamine

Description	MEM Vitamins 100X w/o L-Glutamine
Format	100 ml
Producer	Euroclone
Code Producer	

General Features:

Liquid MEM Vitamins 100X w/o L-Glutamine

pH : 6.5 +/- 1 (before dilution)
Osmolality: as reported in CoA
Endotoxin : < 1 EU/ml
Colour :Yellow

Storage conditions : -20°C in the dark
Shelf life : 24 months

Sterility tests :

- Bacteria in aerobic and anaerobic conditions
- Fungi and yeasts

Concentration of use: 1X

Recommended use :

- Respect storage conditions of the product
- Do not use the product after its expiry date
- Store product in an area protected from light
- Manipulate the product in aseptic conditions (under laminar air flow)
- Wear clothes adapted to the manipulation of the product to avoid contamination

The product is intended to be used in vitro, in laboratory only. Do not use it in therapy, human



or veterinary applications.

Euroclone S.p.A.

A handwritten signature in black ink, appearing to read "Franco Rossetti".