

ECM0749L

Dulbecco's Modified Eagle's Medium Low Glucose without L-Glutamine with Sodium Pyruvate

Description	Dulbecco's Modified Eagle's Medium Low Glucose without L-Glutamine with Sodium Pyruvate
Format	500 ml
Producer	Euroclone
Code Producer	

Specifications:

Liquid medium (1X) suitable for cell cultures containing, among other components, 1000 mg/l di D-Glucose, 110 mg/l Sodium Pyruvate and 3700 mg/l di Sodium Bicarbonate.

Quality control:

Microbiological screening in specific bacteria and mycoplasma mediums.

LAL tested for endotoxins

pH checked with potentiometric measurement at room temperature and osmolality with cryoscopic analysis

Conservation:

4°C for 24 months from production date

Euroclone S.p.A.

