



ECM0920D

Trypsin-EDTA with Phenol Red

Description	Trypsin-EDTA with Phenol Red
Format	100 ml
Producer	Euroclone
Code Producer	

General Features:

The product consists of 0.5 g/l of trypsin porcine and 0.221 g / l of EDTA dissolved in salts without Ca + + and Mg + + ions, with phenol red.

Quality Control

Sterility tests :

- bacteria aerobic-anaerobic condition
- fungi and yeast

pH: 6.8- 7.6

Osmolality : 280 mOsm/kg +/- 10 %

Color: Pink, clear

The enzymatic activity is tested with adherent cells line.

Storage condition:

the solution should be stored at -20 ° C.

Shelf life:

the product under these conditions may be kept for 24 months from date of manufacture.

The product is intended to be used in vitro, in laboratory only. Do not use it in therapy, human or veterinary applications

Euroclone S.p.A.



Thomas Boman



BioNordika Sweden AB
Berzelius Väg 13, 171 65 Solna, Sweden
Tel: +46 (0)8 30 60 10 info@bionordika.se www.bionordika.se