

ECM2001L

RPMI 1640 with stable L-Glutamine

Description	RPMI 1640 with stable L-Glutamine
Format	500 ml
Producer	Euroclone
Code Producer	

General Features:

Liquid medium for cultivation of adherent cells and / or in suspension with stable L-Glutamine

pH : 7.3 ± 0.3

Osmolality: 278 mOsm/kg $\pm 10\%$

Colour: Orange solution

Endotoxin : < 1 EU/ml

Storage conditions : $+2^{\circ}\text{C}$ to $+8^{\circ}\text{C}$ in the dark

Shelf life : 24 months

Sterility tests :

- Bacteria in aerobic and anaerobic conditions
- Fungi and yeasts

Cell growth test :

Medium tested for the ability to support SP2/0-Ag14 cell growth

Recommended use :

- Respect storage conditions of the product
- Do not use the product after its expiry date
- Store product in an area protected from light
- Manipulate the product in aseptic conditions (under laminar air flow)
- Wear clothes adapted to the manipulation of the product to avoid contamination

The product is intended to be used in vitro, in laboratory only. Do not use it in therapy, human or veterinary applications.



Euroclone S.p.A.

Franco Rusconi



BioNordika Sweden AB
Berzelius Väg 13, 171 65 Solna, Sweden
Tel: +46 (0)8 30 60 10 info@bionordika.se www.bionordika.se