



Anti-CAR Linker Antibodies

Detect Surface-Expressed CARs


Make Your CARs GLOW

Challenges with Other CAR Detection Reagents

- Lack of versatility
- Non-specific

Advantages of CST[®] Anti-CAR Linker Antibodies

- Designed to recognize a broad range of scFv-based CARs
- Specifically detect either:
G4S Linker
Whitlow Linker
- Compatible with multiparametric flow panels
- Facilitate the monitoring of CAR expression, trafficking, and persistence
- Lot reservations and bulk orders available



CST offers solutions for more than basic research.

We are eager to help drive your products to market! Let's talk science.



To learn more, please visit:
www.cst-science.com/CARlinkers

 **Cell Signaling**
TECHNOLOGY[®]

Anti-CAR Linker Antibody Reagents

Key CAR Targets

Unconjugated Antibodies

- G4S Linker (E7O2V) Rabbit mAb #71645
- G4S Linker (E7O2V) Rabbit mAb (BSA and Azide Free) #63670
- Whitlow/218 Linker (E3U7Q) Rabbit mAb #57710
- Whitlow/218 Linker (E3U7Q) Rabbit mAb (BSA and Azide Free) #66159

Conjugates

- G4S Linker (E7O2V) Rabbit mAb (PE Conjugate) #38907
- G4S Linker (E7O2V) Rabbit mAb (Alexa Fluor® 647 Conjugate) #69782
- G4S Linker (E7O2V) Rabbit mAb (Alexa Fluor® 488 Conjugate) #50515
- Whitlow/218 Linker (E3U7Q) Rabbit mAb (PE Conjugate) # 62405
- Whitlow/218 Linker (E3U7Q) Rabbit mAb (Alexa Fluor® 647 Conjugate) #69310
- Whitlow/218 Linker ((E3U7Q) Rabbit mAb Alexa Fluor® 488 Conjugate) #55809

PRODUCTS

CD5 (E8X3S) XP® Rabbit mAb #39300

CD7 (E4G1Q) XP® Rabbit mAb #32814

CD19 (Intracellular Domain) (D4V4B) XP® Rabbit mAb #90176

CD20 (E7B7T) XP® Rabbit mAb #48750

CD22 (E7L6Z) XP® Rabbit mAb #98035

CD37 (E4K2M) Rabbit mAb #46894

CD38 (E7Z8C) XP® Rabbit mAb #51000

Claudin-6 (E7U2O) XP® Rabbit mAb #18932

CLEC12A/CLL-1 (E5T5V) Rabbit mAb #57950

CRACC/SLAMF7/CD319 (E5C4M) Rabbit mAb #98611

EGF Receptor vIII (D6T2Q) XP® Rabbit mAb #64952

GUCY2C (E6I2U) Rabbit mAb #37517

HER2/ErbB2 (D8F12) XP® Rabbit mAb #4290

IL-13RA2/CD213a2 (E7U7B) Rabbit mAb #85677

Mesothelin (D9R5G) XP® Rabbit mAb #99966

NG2/CSPG4 (E3B3G) XP® Rabbit mAb #43916

Prostate Specific Membrane Antigen (D7I8E) XP® Rabbit mAb #12815

TACI/TNFRSF13B/CD267 (E1D1P) Rabbit mAb #96641

TNFRSF8/CD30 (E4L4I) XP® Rabbit mAb #54535

TNFRSF8/CD30 (E7E4D) XP® Rabbit mAb #25114

TNFRSF17/BCMA (E6D7B) Rabbit mAb #88183

TRBC1/TCRβ constant region 1 (E6Z3S) Rabbit mAb #79485

ADDITIONAL CAR DETECTION REAGENT

Protein L (PE Conjugate) #58036

