

General Equipment

Below is a selection of instruments offered by BioNordika Sweden AB

If you are looking for anything in particular please email info@bionordika.se

Electrophoresis System

multiSUB Horizontal Gel Systems



- Different options to fit your experiment
- Run up to 210 samples
- Low buffer volumes
- Multichannel pipette compatible combs for speed loading

omniPAGE Vertical Protein Electrophoresis System



- Complete packages for PAGE, Blotting and 2D electrophoresis
- Modular Inserts allow use of a single buffer tank
- Available with Power Supplies

Pipetter



omniPETTE

Choose between:

- Art no: CV10, volume; 0.5 - 10ul
- Art no: CV100, volume; 10 - 100ul
- Art no: CV1000, volume; 100 - 1000ul
- Art no: CV10000, volume; 1000 - 10,000ul
- Art no: CV2, volume; 0.2 - 2ul
- Art no: CV20, volume; 2 - 20ul
- Art no: CV200, volume; 20 - 200ul
- Art no: CV50, volume; 5 - 50ul
- Art no: CV5000, volume; 1000 - 5000ul
- Also available as a pack including 4 and Stand

- Easy operation for right- and left-handed users
- Low pipetting forces
- Highly durable shaft
- UV resistant
- Recalibration system
- Autoclaveable
- Three year warranty



omniPETTE Starter Pack with 2 – 1000 µl pipettes and stand

Art no: CVSTARTERPACK

- Includes CV2, CV20, CV200, CV1000 and CV-4POS Stand
- Easy operation for right- and left-handed users
- Low pipetting forces
- Highly durable shaft
- UV resistant
- Autoclaveable
- Three year warranty



OmniPET Pipette Filler (Pipett boy)

Art no: OMNIPET

- Suitable for 0.5 ml to 100 ml pipettes
- Ergonomically shaped handle
- Sensitive valves for precise work with low volume pipettes
- Protected by filter and safety valve
- Autoclavable nosepiece and pipette holder
- Charging stand optional accessory



Multichannel Pipettes 8 or 12 Channel

omniPETTE 8 Channel or 12 Channel and different volumes; 0,5-10ul, 20-200ul, 5-50ul, 50-300ul

Art no: CV8 or CV12

- 8 Multichannel version
- Easy operation for right- and left-handed users
- Low pipetting forces
- Highly durable shaft
- UV resistant
- Recalibration system
- Autoclaveable
- Three Year Warranty

Vortex



Mini Vortex Mixer

Art no: MINIVORTEX

- ABS and silicone materials to resist most commonly used chemicals
- Low profile design and small footprint (only 4.3" wide) save valuable space on a bench top and in fume hoods
- Rapid vortex action saves time and enables quick mixing of samples
- Modern design allows easy one handed touch start operation and shaking stops when the tube is lifted
- When performing repetitive procedures, a low activation force helps prevent user fatigue
- Maintenance free design and 3 year warranty helps reduce the cost of ownership
- Stable platform reduces the chance of dropping a sample during operation
- Quiet operation and low vibration helps users stay focused

PCR Thermal Cyclers



GTC96S, 96 well Thermal Cycler

Art No: GTC96

- Compatible with 96-well plates, 0.2ml tubes and tube strips
- Protocol optimization selectable from 1 to 24°C across the entire temperature control range from 4-99°C
- Precision temperature control increases both speed and efficiency

Vattenbad



Stirring Water Bath

Available at 10 or 20 liters

Art no: SWB-10L or SWB-20L

- Microprocessor control with digital performance for precise, accurate control
- Built-in magnetic stirring mechanism for water agitation
- Stirring speed control
- LCD screen shows temperature and timer simultaneously
- Safety devices: warning indication on screen to avoid overheating and automatic shut down
- This unit has a 12 month limited warranty