

SPECIAL PRICE

Universal ELISA for Hormones & Small Molecule Detection



Analyte Detection with Universal ELISA Kits

SPECIAL PRICE - UNIVERSAL ELISA FOR HORMONES, CHEMICALS & SMALL MOLECULES

Get a SPECIAL PRICE on any Universal ELISA for prices email info@bionordika.se - use GenieQ124 promo code when placing the order.

VALIDATED

Pre-coated validated ELISA with spike-recovery, linearity, range, CV% & standard curve data.

MULTI-SPECIES

Universal ELISA kits that cross-react with every species including Human, Mouse, NHP, Rat, Bovine, Porcine, & more.

SENSITIVE

Ready-to-use, pre-coated 96-well ELISA with picogram sensitivity for many analytes.

Product Name	Code	Size	Sensitivity	Range	Reactivity	List Price	Promo Price
15d-PGJ2 ELISA	UNES00055	96 assays	1.88ng/mL	3.13-200ng/mL	All Species		
25-HVD3 (25-Hydroxy Vitamin D3) ELISA	UNES00011	96 assays	1.88ng/mL	3.13-200ng/mL	All Species		
3-NT (3-Nitrotyrosine) ELISA	UNES00025	96 assays	0.94ng/mL	1.56-100ng/mL	All Species		
5-HETE (5-Hydroxyeicosatetraenoic Acid) ELISA	UNES00062	96 assays	0.09ng/mL	0.16-10ng/mL	All Species		
5-HIAA (5-Hydroxyindoleacetic Acid) ELISA	UNES00048	96 assays	0.94ng/mL	1.56-100ng/mL	All Species		
6-keto-PGF1a (6-keto-prostaglandin F1a) ELISA	UNES00035	96 assays	9.38pg/mL	15.63-1000pg/mL	All Species		
6-OHDA (6-Hydroxydopamine) ELISA	UNES00065	96 assays	0.09ng/mL	0.16-10ng/mL	All Species		
8-epi-PGF2 alpha (8-Epi-Prostaglandin F2 Alpha) ELISA	UNES00026	96 assays	9.38pg/mL	15.63-1000pg/mL	All Species		
8-OHdG (8-Hydroxydeoxyguanosine) ELISA	UNES00018	96 assays	0.94ng/mL	1.56-100ng/mL	All Species		
Arachidonic Acid ELISA	UNES00032	96 assays	0.94ng/mL	1.56-100ng/mL	All Species		

Product Name	Code	Size	Sensitivity	Range	Reactivity	List Price	Promo Price
Acetylcholine ELISA Kit	UNES00050	96 assays	9.38pg/mL	15.63-1000pg/mL	All Species		
ADMA (Asymmetrical Dimethylarginine) ELISA	UNES00027	96 assays	9.38ng/mL	15.63-1000ng/mL	All Species		
Androsterone ELISA	UNES00047	96 assays	0.47ng/mL	0.78-50ng/mL	All Species		
Aldosterone ELISA	UNES00045	96 assays	18.75pg/mL	31.25-2000pg/mL	All Species		
Bilirubin ELISA Kit	UNES00049	96 assays	0.19Mug/mL	0.31-20Mug/mL	All Species		
Dihydrobiopterin (BH2) ELISA	UNES00060	96 assays	18.75pg/mL	31.25-2000pg/mL	All Species		
cAMP (Cyclic AMP) ELISA	UNES00036	96 assays	0.94ng/mL	1.56-100ng/mL	All Species		
Carnosine ELISA	UNES00044	96 assays	9.38ng/mL	15.63-1000ng/mL	All Species		
cGMP (Cyclic GMP) ELISA	UNES00052	96 assays	0.47pmol/mL	0.78-50pmol/mL	All Species		
Creatinine ELISA	UNES00038	96 assays	0.75Mug/mL	1.25-80Mug/mL	All Species		
Dopamine ELISA	UNES00029	96 assays	18.75pg/mL	31.25-2000pg/mL	All Species		
DHT (Dihydrotestosterone) ELISA	UNES00021	96 assays	23.44pg/mL	39.06-2500pg/mL	All Species		
DHVD3 (1,25-Dihydroxyvitamin D3) ELISA	UNES00012	96 assays	4.69pg/mL	7.81-500pg/mL	All Species		
Epinephrine (Adrenaline) ELISA	UNES00028	96 assays	18.75pg/mL	31.25-2000pg/mL	All Species		
Folic Acid (Vitamin B9) ELISA	UNES00006	96 assays	0.94ng/mL	1.56-100ng/mL	All Species		
GnRH (Gonadotropin-Releasing Hormone) ELISA	UNES00046	96 assays	9.38pg/mL	15.63-1000pg/mL	All Species		
GSH (Glutathione) ELISA	UNES00017	96 assays	0.94Mug/mL	1.56-100Mug/mL	All Species		
Histamine ELISA	UNES00022	96 assays	0.94ng/mL	1.56-100ng/mL	All Species		
Inositol Triphosphate (IP3) ELISA	UNES00039	96 assays	9.38pg/mL	15.63-1000pg/mL	All Species		
Leukotriene B4 (LTB4) ELISA	UNES00041	96 assays	9.38pg/mL	15.63-1000pg/mL	All Species		
Leukotriene E4 (LTE4) ELISA	UNES00061	96 assays	18.75pg/mL	31.25-2000pg/mL	All Species		
Lipoxin A4 (LXA4) ELISA	UNES00034	96 assays	0.47ng/mL	0.78-50ng/mL	All Species		
MDA (Malondialdehyde) ELISA	UNES00040	96 assays	18.75ng/mL	31.25-2000ng/mL	All Species		
Methionine-Enkephalin (MEK) ELISA	UNES00015	96 assays	75pg/mL	125-8000pg/mL	All Species		
Noradrenaline (Norepinephrine) ELISA	UNES00030	96 assays	0.19ng/mL	0.31-20ng/mL	All Species		
Oxytocin ELISA	UNES00019	96 assays	9.38pg/mL	15.63-1000pg/mL	All Species		
Progesterone ELISA Kit	UNES00066	96 assays	0.15ng/mL	0.31-20ng/mL	All Species		
Prostaglandin D2 (PGD2) ELISA	UNES00043	96 assays	9.38pg/mL	15.63-1000pg/mL	All Species		
Prostaglandin E1 (PGE1) ELISA	UNES00033	96 assays	46.88pg/mL	78.13-5000pg/mL	All Species		
Prostaglandin E2 (PGE2) ELISA	UNES00024	96 assays	18.75pg/mL	31.25-2000pg/mL	All Species		
Prostacyclin ELISA	UNES00016	96 assays	18.75pg/mL	31.25-2000pg/mL	All Species		
Pregnenolone ELISA	UNES00054	96 assays	0.47ng/mL	0.78-50ng/mL	All Species		
Pentosidine (PTD) ELISA	UNES00057	96 assays	0.47ng/mL	0.78-50ng/mL	All Species		
Serotonin (5-HT) ELISA	UNES00023	96 assays	9.38ng/mL	15.63-1000ng/mL	All Species		
STAT-3 ELISA Kit	UNES00063	96 assays	0.19ng/mL	0.31-20ng/mL	All Species		
Thromboxane A2 (TXA2) ELISA	UNES00037	96 assays	9.38pg/mL	15.63-1000pg/mL	All Species		
Vitamin B6 ELISA	UNES00005	96 assays	0.94ng/mL	1.56-100ng/mL	All Species		
Vitamin D ELISA	UNES00009	96 assays	3.75ng/mL	6.25-400ng/mL	All Species		
Vitamin D3 ELISA	UNES00010	96 assays	0.94ng/mL	1.56-100ng/mL	All Species		
Vitamin K1 ELISA	UNES00014	96 assays	9.38pg/mL	15.63-1000pg/mL	All Species		